

**SPECTRO® | JET-CARE®**

# webECHO™

[spectro-oil.com](http://spectro-oil.com) | [jet-care.com](http://jet-care.com)



# WELCOME TO webECHO™



Welcome to webECHO™ - the Spectro | Jet-Care secure online portal that gives you immediate access to your latest sample results and trend reports. Designed to work on both desktop and mobile devices, webECHO™ is intuitive and easy to use allowing you to monitor your entire Spectro | Jet-Care analysis and results in one place.

# GET TO KNOW webECHO

## ANALYSIS SUMMARY & DASHBOARD



Analysis Summary



Dashboard

The Analysis Summary provides an overview of the status of your equipment/aircraft fleet enabling you to identify those requiring immediate attention.

Access the dashboard for a particular analysis type by using the tabs at the top of the screen.

The Dashboard provides an overview of the current trend monitoring status of all equipment/aircraft for the selected analysis type.

You can also access the main functions of webECHO as listed:

- 🔍 Operators – view details of all operators.
- 🔍 Equipment/Aircraft – view details of all equipment or aircraft.
- 🔍 Reporting Contacts – view details of all contacts that are defined to receive analysis reports.
- 🔍 Categorise – shows equipment listed by status level e.g. Advanced/Early/Normal/Undefined.
- 🔍 Samples – view details of analysis samples according to analysis type and per equipment.
- 🔍 Find – the operator and equipment icons allow quick filtering access using the drop down sections.
- 🔍 Search – free text search for finding information within the analysis type.

MANAGE, REVIEW  
AND ANALYSE YOUR  
TREND MONITORING  
DATA 24/7/365


# MAIN FUNCTION SCREENS



Operator – View details of all operators



Equipment/Aircraft – view details of all equipment or aircraft.



Categorise – view equipment listed by status level e.g. Advanced/Early/Normal/Undefined.



Laboratory Analysis Graph



Gas Path Analysis (GPA) Graph



Data History – when viewing trend plots you can access the last 50 data points used for analysis for GPA data and all samples for Laboratory data.

# webECHO™

## FUNCTIONALITY

- View and download your PDF reports, data and statistics for:
  - Oil, Hydraulic Fluid, Coolant, Transformer and Fuel Analysis
  - Debris and Filter Analysis by Scanning Electron Microscope
  - Jet-Care Engine Trending by Gas Path Analysis (GPA)
- Use the equipment navigator for quick access to your latest analysis reports, with the ability to graph data for the chosen equipment
- Filter data by alert level, analysis type or any other searchable parameter
- View individual data points in tabular format for quick numerical comparison
- Review the online trends and the analysis of your equipment in real time as well as export the raw data to produce your own graphs and distributions
- Set regular laboratory analysis sample intervals and view 'next sample date' to identify when equipment/aircraft sampling is overdue
- Use the 'Watch' function to keep track of individual equipment or aircraft
- Upload your own notes and documents to individual equipment or aircraft records
- Access diagnostic comments, limits and sample submissions for Laboratory customers
- Access diagnostic comments, data submissions and fault codes for GPA customers
- View detailed images of debris and filter samples taken from the Scanning Electron Microscope
- Use the comprehensive statistics and distributions facility to analyse and interpret laboratory sample data
- Review and request changes to the information held on your organisation
  - locations, operators, equipment, reporting and test suites



# CONTACT

## SPECTRO | JET-CARE



### JET-CARE®

Cedar Knolls | NJ | USA

T +1 973 292 9597

E [webECHO@jet-care.com](mailto:webECHO@jet-care.com)

### SPECTRO® | JET-CARE®

Odiham | Hampshire | UK

T +44 (0) 1256 704000

E [webECHO@spectro-oil.com](mailto:webECHO@spectro-oil.com)

E [webECHO@jet-care.com](mailto:webECHO@jet-care.com)

### SPECTRO OIL AG

Kaiseraugst | Switzerland

T +41 (0) 61 815 90 20

E [webECHO@spectro-oil.com](mailto:webECHO@spectro-oil.com)



[spectro-oil.com](http://spectro-oil.com) | [jet-care.com](http://jet-care.com)